

Product datasheet for **GA204402**

Mouse Tnks activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnks CRISPRa kit - CRISPR gene activation of mouse tankyrase, TRF1-interacting ankyrin-related ADP-ribose polymerase
Symbol:	Tnks
Components:	GA204402G1, Tnks gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204402G2, Tnks gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204402G3, Tnks gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PFX9</u>
Synonyms:	4930554K12Rik; AI662855; ARTD5; C86528; D130072O21Rik; mTNKS1; TANK1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)