

Product datasheet for **GA204403**

Mouse Tnni1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnni1 CRISPRa kit - CRISPR gene activation of mouse troponin I, skeletal, slow 1
Symbol:	Tnni1
UniProt ID:	Q9WUZ5
Synonyms:	2700018B22Rik; AI747285; ssTnl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)