

## Product datasheet for **GA204405**

### Mouse Tnni3 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tnni3 CRISPRa kit - CRISPR gene activation of mouse troponin I, cardiac 3
Symbol:	Tnni3
UniProt ID:	<a href="#">P48787</a>
Synonyms:	cTnI; Tn1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)