

## Product datasheet for **GA204406**

### Mouse Tnnt1 activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tnnt1 CRISPRa kit - CRISPR gene activation of mouse troponin T1, skeletal, slow
Symbol:	Tnnt1
UniProt ID:	<u><a href="#">O88346</a></u>
Synonyms:	AW146156; ssTnT; sTnT; Tnt
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)