

## Product datasheet for **GA204408**

### Mouse Tnnt3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Tnnt3 CRISPRa kit - CRISPR gene activation of mouse troponin T3, skeletal, fast |
| Symbol:       | Tnnt3   |
| UniProt ID:   | <u><a href="#">Q9QZ47</a></u>   |
| Synonyms:     | fTnT  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)         |

#### Product images:



[View online >](#)