

Product datasheet for **GA204410**

Mouse Tnr activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tnr CRISPRa kit - CRISPR gene activation of mouse tenascin R
Symbol:	Tnr
UniProt ID:	Q8BYI9
Synonyms:	J1-tenascin; janusin; restrictin; TN-R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)