

## Product datasheet for **GA204415**

### Mouse Top3a activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Top3a CRISPRa kit - CRISPR gene activation of mouse topoisomerase (DNA) III alpha
Symbol:	Top3a
UniProt ID:	<u><a href="#">O70157</a></u>
Synonyms:	MGC106383
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)