

Product datasheet for **GA204417**

Mouse Tmem165 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem165 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 165
Symbol:	Tmem165
Components:	GA204417G1, Tmem165 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204417G2, Tmem165 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204417G3, Tmem165 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P52875</u>
Synonyms:	AV026557; pFT27; Tpardl; Tparl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)