

Product datasheet for **GA204429**

Mouse Tpst1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tpst1 CRISPRa kit - CRISPR gene activation of mouse protein-tyrosine sulfotransferase 1
Symbol:	Tpst1
Components:	GA204429G1, Tpst1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204429G2, Tpst1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204429G3, Tpst1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O70281</u>
Synonyms:	R75054; Tango13a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)