

Product datasheet for **GA204436**

Mouse Traf3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Traf3 CRISPRa kit - CRISPR gene activation of mouse TNF receptor-associated factor 3
Symbol:	Traf3
UniProt ID:	Q60803
Synonyms:	AI528849; amn; CAP-1; CD40bp; CRAF1; LAP1; T-BAM
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)