

Product datasheet for **GA204445**

Mouse Tfrc activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tfrc CRISPRa kit - CRISPR gene activation of mouse transferrin receptor
Symbol:	Tfrc
UniProt ID:	Q62351
Synonyms:	2610028K12Rik; CD71; E430033M20Rik; Mtvr1; p90; TFR; TFR1; TR; Trfr
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)