

## Product datasheet for **GA204454**

### Mouse *Trpc3* activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	<i>Trpc3</i> CRISPRa kit - CRISPR gene activation of mouse transient receptor potential cation channel, subfamily C, member 3
Symbol:	<i>Trpc3</i>
Synonyms:	Mwk; <i>Trcp3</i> ; <i>Trp3</i> ; <i>Trrp3</i>
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)