

Product datasheet for **GA204475**

Mouse Tsks activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsks CRISPRa kit – CRISPR gene activation of mouse testis-specific serine kinase substrate
Symbol:	Tsks
Components:	GA204475G1, Tsks gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204475G2, Tsks gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204475G3, Tsks gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O54887</u>
Synonyms:	Stk22sl; Tssksl
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

