

Product datasheet for **GA204478**

Mouse Gfus activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tsta3 CRISPRa kit - CRISPR gene activation of mouse tissue specific transplantation antigen P35B
Symbol:	Gfus
UniProt ID:	P23591
Synonyms:	AI256181; FX; Tstap35b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)