

Product datasheet for **GA204479**

Mouse Psm3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psm3 CRISPRa kit - CRISPR gene activation of mouse proteasome (prosome, macropain) 26S subunit, non-ATPase, 3
Symbol:	Psm3
UniProt ID:	P14685
Synonyms:	AI255837; AntP91a; D11Bwg1349e; Psd3; Tstap91a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)