

Product datasheet for **GA204573**

Mouse Vmn2r37 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Vmn2r37 CRISPRa kit - CRISPR gene activation of mouse vomeronasal 2, receptor 37
Symbol:	Vmn2r37
Components:	GA204573G1, Vmn2r37 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204573G2, Vmn2r37 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204573G3, Vmn2r37 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	V2r14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

