

Product datasheet for **GA204637**

Mouse Wnt3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wnt3 CRISPRa kit - CRISPR gene activation of mouse wingless-type MMTV integration site family, member 3
Symbol:	Wnt3
UniProt ID:	P17553
Synonyms:	Int-4; Wnt-3
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)