

Product datasheet for **GA204638**

Mouse Wnt3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wnt3a CRISPRa kit - CRISPR gene activation of mouse wingless-type MMTV integration site family, member 3A
Symbol:	Wnt3a
UniProt ID:	P27467
Synonyms:	vt; Wnt-3a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

