


## Product datasheet for **GA204655**

### Mouse Xlr3c activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Xlr3c CRISPRa kit - CRISPR gene activation of mouse X-linked lymphocyte-regulated 3C
Symbol:	Xlr3c
Components:	GA204655G1, Xlr3c gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204655G2, Xlr3c gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204655G3, Xlr3c gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Xlr3b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)