

Product datasheet for **GA204669**

Mouse Slc23a3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc23a3 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 23 (nucleobase transporters), member 3
Symbol:	Slc23a3
Components:	GA204669G1, Slc23a3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204669G2, Slc23a3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204669G3, Slc23a3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q60850</u>
Synonyms:	SVCT3; YsplI
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

