

Product datasheet for **GA204685**

Mouse Mkrn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mkrn3 CRISPRa kit - CRISPR gene activation of mouse makorin, ring finger protein, 3
Symbol:	Mkrn3
Components:	GA204685G1, Mkrn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204685G2, Mkrn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204685G3, Mkrn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q60764
Synonyms:	D7H15S9-1; Zfp127
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)