

Product datasheet for **GA204710**

Mouse Zfp46 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Zfp46 CRISPRa kit – CRISPR gene activation of mouse zinc finger protein 46
Symbol:	Zfp46
Components:	GA204710G1, Zfp46 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204710G2, Zfp46 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204710G3, Zfp46 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BPP0</u>
Synonyms:	AI426962; Mlz-4; Zfp-46; Znf436
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

