

Product datasheet for **GA204840**

Mouse *Tenm4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tenm4 CRISPRa kit - CRISPR gene activation of mouse teneurin transmembrane protein 4
Symbol:	Tenm4
UniProt ID:	Q3UHK6
Synonyms:	Doc4; ELM2; I(7)-3Rn; I7Rn3; mKIAA1302; Odz4; R75022; Ten-m4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)