

Product datasheet for **GA204906**

Mouse Best1 activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Best1 CRISPRa kit - CRISPR gene activation of mouse bestrophin 1 |
| Symbol: | Best1 |
| UniProt ID: | <u>O88870</u> |
| Synonyms: | Bmd; mBest1; Vmd2 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)