

## Product datasheet for **GA204941**

### Mouse Mapkbp1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Mapkbp1 CRISPRa kit - CRISPR gene activation of mouse mitogen-activated protein kinase binding protein 1                                                                                                                                                                                                         |
| Symbol:       | Mapkbp1                                                                                                                                                                                                                                                                                                          |
| Components:   | GA204941G1, Mapkbp1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA204941G2, Mapkbp1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA204941G3, Mapkbp1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q6NS57</a></u>                                                                                                                                                                                                                                                                                    |
| Synonyms:     | 2810483F24Rik; AW123212; JNKBP-1; Jnkbp1; mKIAA0596                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)