

Product datasheet for **GA204977**

Mouse Psg16 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Psg16 CRISPRa kit – CRISPR gene activation of mouse pregnancy specific glycoprotein 16
Symbol:	Psg16
Components:	GA204977G1, Psg16 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204977G2, Psg16 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204977G3, Psg16 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	bCEA; Ceal1; PSG
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

