

Product datasheet for **GA204985**

Mouse Poli activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Poli CRISPRa kit - CRISPR gene activation of mouse polymerase (DNA directed), iota
Symbol:	Poli
Components:	GA204985G1, Poli gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA204985G2, Poli gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA204985G3, Poli gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6R3M4</u>
Synonyms:	Rad30b
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)