

Product datasheet for **GA205050**

Mouse Polr3e activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Polr3e CRISPRa kit - CRISPR gene activation of mouse polymerase (RNA) III (DNA directed) polypeptide E
Symbol:	Polr3e
Components:	GA205050G1, Polr3e gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205050G2, Polr3e gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205050G3, Polr3e gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9CZT4</u>
Synonyms:	RPC5; Sin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)