

Product datasheet for **GA205054**

Mouse Serinc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Serinc3 CRISPRa kit - CRISPR gene activation of mouse serine incorporator 3
Symbol:	Serinc3
Components:	GA205054G1, Serinc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205054G2, Serinc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205054G3, Serinc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QZL9</u>
Synonyms:	AA959945; AA960047; AIGP1; DIFF33; Tde1; TMS-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)