

## Product datasheet for **GA205091**

### Mouse Eif2ak4 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Eif2ak4 CRISPRa kit - CRISPR gene activation of mouse eukaryotic translation initiation factor 2 alpha kinase 4 |
| Symbol:       | Eif2ak4   |
| UniProt ID:   | <a href="#">Q9QZ05</a>  |
| Synonyms:     | 2610011M03; GCN2; MGCN2   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)   |

#### Product images:



[View online >](#)