

Product datasheet for **GA205092**

Mouse Tlx3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tlx3 CRISPRa kit - CRISPR gene activation of mouse T cell leukemia, homeobox 3
Symbol:	Tlx3
Components:	GA205092G1, Tlx3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205092G2, Tlx3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205092G3, Tlx3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>O55144</u>
Synonyms:	Hox11l2; Rnx; Tlx1l2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)