


Product datasheet for **GA205107**

Mouse Ddx24 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ddx24 CRISPRa kit - CRISPR gene activation of mouse DEAD (Asp-Glu-Ala-Asp) box polypeptide 24
Symbol:	Ddx24
Components:	GA205107G1, Ddx24 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205107G2, Ddx24 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205107G3, Ddx24 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9ESV0</u>
Synonyms:	1700055J08Rik; 2510027P10Rik; AI649272
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)