

Product datasheet for **GA205108**

Mouse Pla2g7 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pla2g7 CRISPRa kit - CRISPR gene activation of mouse phospholipase A2, group VII (platelet-activating factor acetylhydrolase, plasma)
Symbol:	Pla2g7
Components:	GA205108G1, Pla2g7 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205108G2, Pla2g7 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205108G3, Pla2g7 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q60963</u>
Synonyms:	PAF-AH; R75400
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)