

Product datasheet for **GA205114**

Mouse Golga5 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Golga5 CRISPRa kit - CRISPR gene activation of mouse golgi autoantigen, golgin subfamily a, 5
Symbol:	Golga5
Components:	GA205114G1, Golga5 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205114G2, Golga5 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205114G3, Golga5 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QYE6</u>
Synonyms:	Ret-II
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)