

Product datasheet for **GA205116**

Mouse Phlda3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Phlda3 CRISPRa kit - CRISPR gene activation of mouse pleckstrin homology like domain, family A, member 3
Symbol:	Phlda3
UniProt ID:	Q9WV95
Synonyms:	Tih1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)