

Product datasheet for **GA205153**

Mouse Skp2 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | Skp2 CRISPRa kit - CRISPR gene activation of mouse S-phase kinase-associated protein 2 (p45) |
| Symbol: | Skp2 |
| UniProt ID: | Q9Z0Z3 |
| Synonyms: | 4930500A04Rik; FBXL1; FWD1; p45 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)