

Product datasheet for **GA205165**

Mouse Sergef activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sergef CRISPRa kit - CRISPR gene activation of mouse secretion regulating guanine nucleotide exchange factor
Symbol:	Sergef
Components:	GA205165G1, Sergef gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205165G2, Sergef gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205165G3, Sergef gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80YD6</u>
Synonyms:	DelGEF; Gef; Gnef; Gnefr
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)