

Product datasheet for **GA205183**

Mouse Tada1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tada1 CRISPRa kit - CRISPR gene activation of mouse transcriptional adaptor 1
Symbol:	Tada1
UniProt ID:	Q99LM9
Synonyms:	2900026B15Rik; D1Ertd251e; Staf42; Tada1l
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)