

Product datasheet for **GA205184**

Mouse Tango2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tango2 CRISPRa kit - CRISPR gene activation of mouse transport and golgi organization 2
Symbol:	Tango2
Components:	GA205184G1, Tango2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205184G2, Tango2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205184G3, Tango2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	D16H22S680E; T10
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)