

Product datasheet for **GA205189**

Mouse Imp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Imp4 CRISPRa kit - CRISPR gene activation of mouse IMP4, U3 small nucleolar ribonucleoprotein
Symbol:	Imp4
Components:	GA205189G1, Imp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205189G2, Imp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205189G3, Imp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8VHZ7</u>
Synonyms:	AA409888; AV031295; D1Wsu40e
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)