

Product datasheet for **GA205198**

Mouse Selenoi activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | Selenoi CRISPRa kit - CRISPR gene activation of mouse selenoprotein I |
| Symbol: | Selenoi |
| UniProt ID: | <u>Q80TA1</u> |
| Synonyms: | 4933402G07Rik; AI448296; AI452230; C79563; D5Wsu178e; Ept1; mKIAA1724; SELI |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)