

Product datasheet for **GA205204**

Mouse Trim36 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim36 CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 36
Symbol:	Trim36
Components:	GA205204G1, Trim36 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205204G2, Trim36 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205204G3, Trim36 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80WG7</u>
Synonyms:	D18Wsu100e; haprin
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)