

Product datasheet for **GA205250**

Mouse *Iqgap1* activation kit by CRISPRa

Product data:

| | |
|---------------|---|
| Product Type: | CRISPRa Kits |
| Description: | <i>Iqgap1</i> CRISPRa kit - CRISPR gene activation of mouse IQ motif containing GTPase activating protein 1 |
| Symbol: | <i>Iqgap1</i> |
| UniProt ID: | Q9JKF1 |
| Synonyms: | AA682088; D7Ertd237e; D7Ertd257e; mKIAA0051 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)