

Product datasheet for **GA205258**

Mouse Pcsk1n activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcsk1n CRISPRa kit - CRISPR gene activation of mouse proprotein convertase subtilisin/kexin type 1 inhibitor
Symbol:	Pcsk1n
Components:	GA205258G1, Pcsk1n gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205258G2, Pcsk1n gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205258G3, Pcsk1n gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q9QXV0</u>
Synonyms:	AI848336; bLEN; KEP; ILEN; Pan3; PEN; PEN19; PEN20; proSAAS; SAAS; SAASCT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)