

Product datasheet for **GA205290**

Mouse Tor3a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tor3a CRISPRa kit - CRISPR gene activation of mouse torsin family 3, member A
Symbol:	Tor3a
UniProt ID:	<u>Q9ER38</u>
Synonyms:	Adir
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)