

Product datasheet for **GA205334**

Mouse Hlf4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hist1h1e CRISPRa kit - CRISPR gene activation of mouse histone cluster 1, H1e
Symbol:	Hlf4
Components:	GA205334G1, Hlf4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205334G2, Hlf4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205334G3, Hlf4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>P43274</u>
Synonyms:	H1-4; H1.4; H1e; Hlf4; H1s-4; H1var2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

