


Product datasheet for **GA205374**

Mouse Sh3bgr activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sh3bgr CRISPRa kit - CRISPR gene activation of mouse SH3-binding domain glutamic acid-rich protein
Symbol:	Sh3bgr
Components:	GA205374G1, Sh3bgr gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA205374G2, Sh3bgr gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA205374G3, Sh3bgr gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9WUZZ
Synonyms:	5430437A18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)