

Product datasheet for **GA205413**

Mouse H2az1 activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | H2azf CRISPRa kit - CRISPR gene activation of mouse H2A histone family, member Z |
| Symbol: | H2az1 |
| UniProt ID: | <u>P0C0S6</u> |
| Synonyms: | H2A.Z; H2a.z-1 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)