

Product datasheet for **GA205418**

Mouse *Srrm1* activation kit by CRISPRa

Product data:

| | |
|---------------|--|
| Product Type: | CRISPRa Kits |
| Description: | <i>Srrm1</i> CRISPRa kit - CRISPR gene activation of mouse serine/arginine repetitive matrix 1 |
| Symbol: | <i>Srrm1</i> |
| UniProt ID: | Q52K18 |
| Synonyms: | AA407769; POP101; <i>Srm160</i> |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)