

## Product datasheet for **GA205433**

### Mouse Knstrn activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Knstrn CRISPRa kit - CRISPR gene activation of mouse kinetochore-localized astrin/SPAG5 binding
Symbol:	Knstrn
UniProt ID:	<a href="#">Q9D9Z1</a>
Synonyms:	1700025D04Rik; C15orf23; D2Ertd750e; SKAP; Traf4af1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online >](#)